

TUBERCULOSIS DETECTION IN CHEST RADIOGRAPHY



MINI PROJECT REPORT

Submitted by

MOHANASUNDAR S [20AD033]

TAMILARASU A [20AD051]

AYYAPPAN M [20AD007]

In partial fulfillment for the award of the degree of

BACHELOR OF TECHNOLOGY

IN

ARTIFICAL INTELLIGENCE AND DATA SCIENCE

MUTHAYAMMAL ENGINEERING COLLEGE(AUTONOMOUS) RASIPURAM

DEC 2022

CHAPTER 6

PROJECT DESCRIPTION

6.1 System Implementation




Fig no: 6.1 Proposed image processing pipeline.

Possible workflow for a teleradiological TB detection service for health care providers working remotely (schematic). After taking a photograph of the original chest X-ray with a smartphone application, the image is transferred wirelessly to a remote server, where image analysis is carried out in two steps: anomaly detection and classification. The server then sends the image with an overlaid heatmap and a structured report back to the user (health care provider), who may use it to augment his/her clinical judgement for optimal management of the patient.




Fig .no:6.2.1 Example of X-ray images and corresponding ground truth lung masks from Kaggle dataset

During the phase of collecting the photographs for the dataset, images with smaller resolution and dimension but 500 pixels weren't considered as valid imagesfor the dataset. additionally, only the pictures where the region of interest was in higher resolution were marked as eligible candidates for the dataset. in this way, it had been ensured that images contain all the needed information for feature learning. Many resources may be found by searching across the net, but their relevance is usually unreliable, within the interest of confirming the accuracy of classes within the dataset, initially grouped by a keywords search, medical experts examined chest x-ray images and labeled all the photographs with appropriate disease acronym.

Because it is thought, it's important to use accurately classified images for thetraining and validation dataset. Only in this way may an appropriate and reliable detecting model be developed. during this stage, duplicated images that were left after the initial iteration of gathering and grouping images into




Fig no 7.1: Architecture of A) original U-Net and B) modified U-Net.

7.2 TB CLASSIFICATION

As mentioned earlier, there are two different experiments (using non-segmented and segmented lungs X-ray images) were conducted for the classification of TB and normal (non-TB) cases. The comparative performance for different CNNs for the binary classification is shown in Table 5. It is apparent from Table 5 that all the evaluated pre-trained models perform very well in classifying TB and normal images in this two-class problem. Among the networks trained with X-ray images without segmentation, CheXNet is performing better for classifying the X-ray images. Even though CheXNet is shallower than DenseNet201, it was originally trained on X-ray




Fig no .7.2.1:Comparison of the ROC curves for Normal, and Tuberculosis classification using CNN based models for non-segmented (A) and segmented(B) CXR images.

Classification Report:							
	precision	recall	f1-score	support			
Normal	0.9917	1.0000	0.9959	120			
Tuberculosis	1.0000	0.9875	0.9937	80			
accupacy			0.9950	200			
accuracy							
macro avg	0.9959	0.9938	0.9948	200			
weighted avg	0.9950	0.9950	0.9950	200			




Fig no: 7.2.2 Confusion matrix for Normal and Tuberculosis (TB) classification for ChexNet model without segmented X-ray (A), and DenseNet201 model with segmented X-ray (B).

SUPPLEMENTARY MATERIALS

A chest X-ray database of 796 Tuberculosis patients' image and 1200 normal images were released. This database was created from the 4 publicly available databases, which are referenced in the database. https://www.kaggle.com/ tawsifurrahman/tuberculosis-tb-chest-x ray-dataset




Fig no.8: Overview of the complete system using methodology in tuberculosis detection system

APPENDIX A SOURCE CODE

```
import pandas as pd
import numpy as np
os.environ['TF CPP MIN LOG LEVEL'] = '2'
import time
import matplotlib.pyplot as plt
import seaborn as sns
sns.set style('darkgrid')
import shutil
from sklearn.metrics import confusion matrix, classification report
from sklearn.model_selection import train_test_split
import tensorflow as tf
from tensorflow import keras
from tensorflow.keras.preprocessing.image import ImageDataGenerator
from tensorflow.keras.layers import Dense, Activation, Dropout, Conv2D
, MaxPooling2D, BatchNormalization
from tensorflow.keras.optimizers import Adam, Adamax
from tensorflow.keras.metrics import categorical crossentropy
from tensorflow.keras import regularizers
from tensorflow.keras.models import Model
from tensorflow.keras import backend as K
import time
from sklearn.metrics import f1 score
from IPython.display import YouTubeVideo
import sys
if not sys.warnoptions:
   warnings.simplefilter("ignore")
pd.set option('display.max columns', None) # or 1000
pd.set option('display.max rows', None) # or 1000
pd.set option('display.max colwidth', None) # or 199
print ('Modules loaded')
def print in color(txt msg, fore tupple=(0,255,255), back tupple=(100,
100,100)):
   rf,gf,bf=fore tupple
    rb, qb, bb=back tupple
   msg='{0}' + txt msg
    mat='\33[38;2;' + str(rf) +';' + str(gf) + ';' + str(bf) + ';48;
2;' + str(rb) + ';' +str(gb) + ';' + str(bb) + 'm'
```

```
print(msg .format(mat), flush=True)
    print('\33[0m', flush=True) # returns default print color to bac
msg='test of default colors'
print in color(msg)
def make dataframes(sdir):
    filepaths=[]
    labels=[]
    classlist=sorted(os.listdir(sdir) )
    for klass in classlist:
        classpath=os.path.join(sdir, klass)
        if os.path.isdir(classpath):
            flist=sorted(os.listdir(classpath))
            desc=f'{klass:25s}'
            for f in tqdm(flist, ncols=130,desc=desc, unit='files',
colour='blue'):
                fpath=os.path.join(classpath,f)
                filepaths.append(fpath)
                labels.append(klass)
    Fseries=pd.Series(filepaths, name='filepaths')
    Lseries=pd.Series(labels, name='labels')
    df=pd.concat([Fseries, Lseries], axis=1)
    train df, dummy df=train test split(df, train size=.8, shuffle=T
rue, random state=123, stratify=df['labels'])
    valid df, test df=train test split(dummy df, train size=.5, shuf
fle=True, random state=123, stratify=dummy df['labels'])
    classes=sorted(train df['labels'].unique())
    class count=len(classes)
    sample df=train df.sample(n=50, replace=False)
    wt=0
    for i in range(len(sample df)):
        fpath=sample df['filepaths'].iloc[i]
            img=cv2.imread(fpath)
            h=img.shape[0]
            w=imq.shape[1]
            wt +=w
            ht +=h
            count +=1
    have=int(ht/count)
```

```
wave=int(wt/count)
    aspect ratio=have/wave
    print('number of classes in processed dataset= ', class count)
    counts=list(train df['labels'].value counts())
    print(counts[0], type(counts[0]))
   print('the maximum files in any class in train df is ', max(coun
   print('train df length: ', len(train df), ' test df length: ',
len(test df), ' valid df length: ', len(valid df))
   print('average image height= ', have, ' average image width= ',
wave, ' aspect ratio h/w= ', aspect ratio)
    return train df, test df, valid df, classes, class count
sdir=r'../content/drive/MyDrive/TB Chest Radiography Database'
train df, test df, valid df, classes, class count=make dataframes(sd
ir)
def trim(df, max samples, min samples, column):
   df=df.copy()
   classes=df[column].unique()
   class count=len(classes)
    length=len(df)
   print ('dataframe initially is of length ',length, ' with ', cla
ss count, ' classes')
    groups=df.groupby(column)
    trimmed df = pd.DataFrame(columns = df.columns)
    groups=df.groupby(column)
    for label in df[column].unique():
        group=groups.get group(label)
        count=len(group)
        if count > max samples:
            sampled group=group.sample(n=max samples, random state=1
23, axis=0)
            trimmed df=pd.concat([trimmed df, sampled group], axis=0
            if count>=min samples:
                sampled group=group
                trimmed df=pd.concat([trimmed df, sampled group], ax
is=0)
   print('after trimming, the maximum samples in any class is now '
,max samples, ' and the minimum samples in any class is ', min sampl
    classes=trimmed df[column].unique()# return this in case some cl
```

```
length=len(trimmed df)
    print ('the trimmed dataframe now is of length ', length, ' with
', class count, ' classes')
    return trimmed df, classes, class count
max samples=500
min samples=500
column='labels'
train df, classes, class count=trim(train df, max samples, min sampl
def balance(df, n, working dir, img size):
    df=df.copy()
    print('Initial length of dataframe is ', len(df))
    aug dir=os.path.join(working dir, 'aug')# directory to store aug
    if os.path.isdir(aug dir):# start with an empty directory
        shutil.rmtree(aug dir)
    os.mkdir(aug dir)
    for label in df['labels'].unique():
        dir path=os.path.join(aug dir,label)
        os.mkdir(dir path) # make class directories within aug direc
    total=0
    gen=ImageDataGenerator(horizontal flip=True, rotation range=20,
 width shift range=.2,
                                  height shift range=.2, zoom range=
    groups=df.groupby('labels') # group by class
    for label in df['labels'].unique(): # for every class
        group=groups.get group(label) # a dataframe holding only ro
        sample count=len(group) # determine how many samples there
        if sample count < n: # if the class has less than target numb
            aug img count=0
           delta=n - sample count # number of augmented images to
           target_dir=os.path.join(aug dir, label) # define where
            msg='{0:40s} for class {1:^30s} creating {2:^5s} augment
ed images'.format(' ', label, str(delta))
            print(msg, '\r', end='') # prints over on the same line
            aug gen=gen.flow from dataframe( group, x col='filepath
s', y col=None, target size=img size,
```

```
class mode=None, batch s
ize=1, shuffle=False,
                                             save to dir=target dir,
save prefix='aug-', color mode='rgb',
                                             save format='jpg')
            while aug img count<delta:</pre>
                images=next(aug gen)
                aug img count += len(images)
            total +=aug img count
    print('Total Augmented images created= ', total)
    aug fpaths=[]
    aug labels=[]
    classlist=os.listdir(aug dir)
    for klass in classlist:
        classpath=os.path.join(aug dir, klass)
        flist=os.listdir(classpath)
        for f in flist:
            fpath=os.path.join(classpath,f)
            aug fpaths.append(fpath)
            aug labels.append(klass)
    Fseries=pd.Series(aug fpaths, name='filepaths')
    Lseries=pd.Series(aug labels, name='labels')
    aug df=pd.concat([Fseries, Lseries], axis=1)
    df=pd.concat([df,aug df], axis=0).reset index(drop=True)
    print('Length of augmented dataframe is now ', len(df))
    return df
def make gens (batch size, train df, test df, valid df, img size):
    trgen=ImageDataGenerator()
    t and v gen=ImageDataGenerator()
    msg='{0:70s} for train generator'.format(' ')
    print(msg, '\r', end='') # prints over on the same line
    train gen=trgen.flow from dataframe(train df, x col='filepaths',
y col='labels', target size=img size,
                                         class mode='categorical', co
lor mode='rgb', shuffle=True, batch size=batch size)
    msg='{0:70s} for valid generator'.format(' ')
    print(msg, '\r', end='') # prints over on the same line
    valid gen=t and v gen.flow from dataframe(valid df, x col='filep
aths', y col='labels', target size=img size,
                                        class mode='categorical', col
or mode='rgb', shuffle=False, batch size=batch size)
```

```
length=len(test df)
    test batch size=sorted([int(length/n) for n in range(1,length+1)
 if length % n ==0 and length/n<=80], reverse=True) [0]</pre>
    test steps=int(length/test batch size)
    msg='{0:70s} for test generator'.format(' ')
    print(msg, '\r', end='') # prints over on the same line
    test gen=t and v gen.flow from dataframe(test df, x col='filepat
hs', y col='labels', target size=img size,
                                       class mode='categorical', col
or mode='rgb', shuffle=False, batch size=test batch size)
    classes=list(train gen.class indices.keys())
    class indices=list(train gen.class indices.values())
    class count=len(classes)
    labels=test gen.labels
   print ( 'test batch size: ' ,test batch size, ' test steps: ',
test steps, ' number of classes : ', class count)
    return train gen, test gen, valid gen, test steps
img size=(224,224)
batch size =30
train gen, test gen, valid gen, test steps =make gens(batch size, tr
ain df, test df, valid df, img size)
def show image samples(gen ):
    t dict=gen.class indices
    classes=list(t dict.keys())
    images,labels=next(gen) # get a sample batch from the generator
    plt.figure(figsize=(25, 25))
    length=len(labels)
    if length<25: #show maximum of 25 images
        r=length
    else:
    for i in range(r):
        plt.subplot(5, 5, i + 1)
        image=images[i] /255
        plt.imshow(image)
        index=np.argmax(labels[i])
        class name=classes[index]
        plt.title(class name, color='blue', fontsize=18)
        plt.axis('off')
    plt.show()
```

```
show image samples(train gen )
def make model(img size, lr, mod num=3):
    img shape=(img size[0], img size[1], 3)
        base model=tf.keras.applications.efficientnet.EfficientNetB0
(include top=False, weights="imagenet",input_shape=img_shape, poolin
q='max')
        msg='Created EfficientNet B0 model'
    elif mod num == 3:
        base model=tf.keras.applications.efficientnet.EfficientNetB3
(include top=False, weights="imagenet", input shape=img shape, poolin
q='max')
        msq='Created EfficientNet B3 model'
    elif mod num == 5:
        base model=tf.keras.applications.efficientnet.EfficientNetB5
(include top=False, weights="imagenet", input shape=img shape, poolin
g='max')
        msg='Created EfficientNet B5 model'
        base model=tf.keras.applications.efficientnet.EfficientNetB7
(include top=False, weights="imagenet", input shape=img shape, poolin
q='max')
        msq='Created EfficientNet B7 model'
    base model.trainable=True
    x=base model.output
    x=BatchNormalization(axis=-1, momentum=0.99, epsilon=0.001)(x)
    x = Dense(256, kernel regularizer = regularizers.12(l = 0.016),a
ctivity regularizer=regularizers.11(0.006),
                    bias regularizer=regularizers.11(0.006) ,activat
ion='relu')(x)
    x=Dropout(rate=.4, seed=123)(x)
    output=Dense(class count, activation='softmax')(x)
    model=Model(inputs=base model.input, outputs=output)
    model.compile(Adamax(learning rate=lr), loss='categorical crosse
ntropy', metrics=['accuracy'])
    msg=msg + f' with initial learning rate set to {lr}'
    print in color(msg)
    return model
lr=.001
model=make model(img size, lr) # using B3 model by default
lass LR ASK(keras.callbacks.Callback):
    def __init__ (self, model, epochs, ask_epoch, dwell=True, facto
```

```
super(LR ASK, self). init ()
        self.model=model
        self.ask epoch=ask epoch
        self.epochs=epochs
        self.ask=True # if True query the user on a specified epoch
        self.lowest vloss=np.inf
        self.lowest aloss=np.inf
        self.best weights=self.model.get weights() # set best weight
       self.best epoch=1
       self.plist=[]
       self.alist=[]
       self.dwell= dwell
        self.factor=factor
   def get list(self): # define a function to return the list of %
        return self.plist, self.alist
   def on train begin(self, logs=None): # this runs on the beginnin
        if self.ask epoch == 0:
           print('you set ask epoch = 0, ask epoch will be set to 1
', flush=True)
            self.ask epoch=1
       if self.ask epoch >= self.epochs: # you are running for epoc
           print('ask epoch >= epochs, will train for ', epochs, '
epochs', flush=True)
            self.ask=False # do not query the user
       if self.epochs == 1:
           msg =f'Training will proceed until epoch {ask epoch} the
           print in color(msg )
           msg='enter H to halt training or enter an integer for ho
           print in color(msg)
            if self.dwell:
               msg='learning rate will be automatically adjusted du
ring training'
                print in color(msg, (0,255,0))
```

```
def on train end(self, logs=None): # runs at the end of traini
        msg=f'loading model with weights from epoch {self.best epoch
        print in color(msg, (0,255,255))
        self.model.set weights(self.best weights) # set the weights
        tr duration=time.time() - self.start time  # determine how
       hours = tr duration // 3600
       minutes = (tr duration - (hours * 3600)) // 60
        seconds = tr duration - ((hours * 3600) + (minutes * 60))
        msg = f'training elapsed time was {str(hours)} hours, {minut
es:4.1f} minutes, {seconds:4.2f} seconds)'
       print in color (msg) # print out training duration time
    def on epoch end(self, epoch, logs=None): # method runs on the
        vloss=logs.get('val loss') # get the validation loss for th
        aloss=logs.get('loss')
        if epoch >0:
            deltav = self.lowest vloss- vloss
            pimprov=(deltav/self.lowest vloss) * 100
           self.plist.append(pimprov)
           deltaa=self.lowest aloss-aloss
           aimprov=(deltaa/self.lowest aloss) * 100
            self.alist.append(aimprov)
           pimprov=0.0
           aimprov=0.0
        if vloss< self.lowest vloss:</pre>
            self.lowest vloss=vloss
            self.best weights=self.model.get weights() # set best we
            self.best epoch=epoch + 1
           msg=f'\n validation loss of {vloss:7.4f} is {pimprov:7.4
f} % below lowest loss, saving weights from epoch {str(epoch + 1):3s
           print in color(msg, (0,255,0)) # green foreground
           pimprov=abs(pimprov)
```

```
msg=f'\n validation loss of {vloss:7.4f} is {pimprov:7.4
f} % above lowest loss of {self.lowest vloss:7.4f} keeping weights f
rom epoch {str(self.best epoch)} as best weights'
            print in color(msg, (255,255,0)) # yellow foreground
            if self.dwell: # if dwell is True when the validation lo
                lr=float(tf.keras.backend.get value(self.model.optim
izer.lr)) # get the current learning rate
                new lr=lr * self.factor
                msg=f'learning rate was automatically adjusted from
{lr:8.6f} to {new lr:8.6f}, model weights set to best weights'
                print in color(msg) # cyan foreground
                tf.keras.backend.set value(self.model.optimizer.lr,
new lr) # set the learning rate in the optimizer
                self.model.set weights(self.best weights) # set the
        if aloss< self.lowest aloss:</pre>
            self.lowest aloss=aloss
        if self.ask: # are the conditions right to query the user?
            if epoch + 1 ==self.ask epoch: # is this epoch the one f
                msg='press enter to continue or enter a comment bel
                print in color(msg)
                comment=input(' ')
                if comment !='':
                    print in color(comment, (155,245,66))
                msg='\n Enter H to end training or an integer for t
                print in color(msg) # cyan foreground
                ans=input()
                    msg=f'you entered {ans}, Training halted on epo
ch {epoch+1} due to user input\n'
                    print in color(msg)
                    self.model.stop training = True # halt training
                    self.ask epoch += int(ans)
                    if self.ask epoch > self.epochs:
```


```
print('\nYou specified maximum epochs of as
', self.epochs, ' cannot train for ', self.ask epoch, flush =True)
                        msg=f'you entered {ans} Training will contin
ue to epoch {self.ask epoch}'
                        print in color(msg) # cyan foreground
                        if self.dwell==False:
                            lr=float(tf.keras.backend.get value(self
.model.optimizer.lr)) # get the current learning rate
                            msg=f'current LR is {lr:8.6f} hit ente
                            print in color(msg) # cyan foreground
                            ans=input(' ')
                            if ans =='':
                                msg=f'keeping current LR of {lr:7.5f
                                print in color(msg) # cyan foregroun
                                new lr=float(ans)
                                tf.keras.backend.set value(self.mode
1.optimizer.lr, new lr) # set the learning rate in the optimizer
                                msg=f' changing LR to {ans}'
                                print in color(msg) # cyan foregroun
epochs=40
ask epoch=12
ask=LR ASK(model, epochs, ask epoch)
callbacks=[ask]
history=model.fit(x=train gen, epochs=epochs, verbose=1, callbacks=
callbacks, validation data=valid gen,
                 shuffle=False, initial epoch=0)
def tr plot(tr data, start epoch):
   tacc=tr data.history['accuracy']
   tloss=tr data.history['loss']
   vacc=tr data.history['val accuracy']
   vloss=tr data.history['val loss']
   Epoch count=len(tacc) + start epoch
   Epochs=[]
    for i in range (start epoch , Epoch count):
        Epochs.append(i+1)
    index loss=np.argmin(vloss)# this is the epoch with the lowest
```


```
val lowest=vloss[index loss]
   index acc=np.argmax(vacc)
   acc highest=vacc[index acc]
   plt.style.use('fivethirtyeight')
   sc label='best epoch= '+ str(index loss+1 +start epoch)
   vc label='best epoch= '+ str(index acc + 1+ start epoch)
   fig, axes=plt.subplots(nrows=1, ncols=2, figsize=(25,10))
   axes[0].plot(Epochs,tloss, 'r', label='Training loss')
   axes[0].plot(Epochs, vloss, 'g', label='Validation loss')
   axes[0].scatter(index loss+1 +start epoch, val lowest, s=150, c=
'blue', label=sc label)
   axes[0].scatter(Epochs, tloss, s=100, c='red')
   axes[0].set title('Training and Validation Loss')
   axes[0].set xlabel('Epochs', fontsize=18)
   axes[0].set ylabel('Loss', fontsize=18)
   axes[0].legend()
   axes[1].plot (Epochs, tacc, 'r', label= 'Training Accuracy')
   axes[1].scatter(Epochs, tacc, s=100, c='red')
   axes[1].plot (Epochs, vacc, 'g', label= 'Validation Accuracy')
   axes[1].scatter(index acc+1 +start epoch,acc_highest, s=150, c=
'blue', label=vc label)
   axes[1].set title('Training and Validation Accuracy')
   axes[1].set xlabel('Epochs', fontsize=18)
   axes[1].set ylabel('Accuracy', fontsize=18)
   axes[1].legend()
   plt.tight layout
   plt.show()
```


APPENDIX B


Screen Shots

Image table in server where the images are uploaded










U-Net output segmentation results for two sample lung CXR images. a shows the original images. b shows the results of U-Net. c show the results after ROI extraction


Clustering Of The Input Image

Matlab Command Prompt Where The Affected Area Is Displayed



Classification Report:							
	precision	recall	f1-score	support			
Normal	0.9917	1.0000	0.9959	120			
Tuberculosis	1.0000	0.9875	0.9937	80			
accuracy			0.9950	200			
macro avg	0.9959	0.9938	0.9948	200			
weighted avg	0.9950	0.9950	0.9950	200			

Help Dialog Displaying the Tuberculosis disease name






Fig no.9: Overview of the complete system in detection system